

# OHAUS 終売品と後継機/同等スペック機

※下記記載の後継機種・同等スペック製品は必ずしも終売品のスペックと完全合致するものではありません。

※下記に記載がない製品についてはオーハスまでお問い合わせをお願いいたします。

更新日:2023年10月23日

終売品シリーズ	型式	ひょう量	最小表示	販売終了	後継機種・同等スペック機型式	品番	ひょう量	最小表示
ボイジャーシリーズ	V10640	62g	0.1mg	2003	EX124G	83021448	120g	0.1mg
	V11140	110g	0.1mg		EX124G	83021448	120g	0.1mg
	V12140	210g	0.1mg		EX224G	83021449	220g	0.1mg
	V1RR80	100g/210g	0.1mg/1mg		EX224G	83021449	220g	0.1mg
	V12130	210g	1mg		EX223G	83021452	220g	0.001g
	V14130	410g	1mg		EX423G	83021453	420g	0.001g
	V1RV70	100g/410g	1mg/0.01g		EX423G	83021453	420g	0.001g
	V16120	610g	0.01g		EX2202G	83021459	2200g	0.01g
	V1B120	2100g	0.01g		EX2202G	83021459	2200g	0.01g
	V1D120	4100g	0.01g		EX4202G	83021460	4200g	0.01g
	V1RW60	1000g/4100g	0.01g/0.1g		EX4202G	83021460	4200g	0.01g
	V1D110	4100g	0.1g		EX6201G	83021462	6200g	0.1g
	V1F110	6100g	0.1g		EX6201G	83021462	6200g	0.1g
	V1H110	8100g	0.1g		EX10201G	83021464	10200g	0.1g
	V00640	62g	0.1mg		EX124G	83021448	120g	0.1mg
	V01140	110g	0.1mg		EX124G	83021448	120g	0.1mg
	V02140	210g	0.1mg		EX224G	83021449	220g	0.1mg
	V0RR80	100g/210g	0.1mg/1mg		EX224G	83021449	220g	0.1mg
	V02130	210g	1mg		EX223G	83021452	220g	0.001g
	V04130	410g	1mg		EX423G	83021453	420g	0.001g
	V0RV70	100g/410g	1mg/0.01g		EX423G	83021453	420g	0.001g
	V06120	610g	0.01g		EX2202G	83021459	2200g	0.01g
	V0B120	2100g	0.01g		EX2202G	83021459	2200g	0.01g
	V0D120	4100g	0.01g		EX4202G	83021460	4200g	0.01g
	V0RW60	1000g/4100g	0.01g/0.1g		EX4202G	83021460	4200g	0.01g
	V0D110	4100g	0.1g		EX6201G	83021462	6200g	0.1g
V0F110	6100g	0.1g	EX6201G	83021462	6200g	0.1g		
V0H110	8100g	0.1g	EX10201G	83021464	10200g	0.1g		
エクスプローラープロシリーズ	EP64C	62g	0.1mg	2010	EX124G	83021448	120g	0.1mg
	EP114C	110g	0.1mg		EX124G	83021448	120g	0.1mg
	EP214C	210g	0.1mg		EX224G	83021449	220g	0.1mg
	EP213C	210g	1mg		EX223G	83021452	220g	0.001g
	EP213	210g	1mg		EX223G	83021452	220g	0.001g
	EP413C	410g	1mg		EX423G	83021453	420g	0.001g
	EP413	410g	1mg		EX423G	83021453	420g	0.001g
	EP613C	610g	1mg		EX623G	83021455	620g	0.001g
	EP613	610g	1mg		EX623G	83021455	620g	0.001g
	EP2102C	2100g	0.01g		EX2202G	83021459	2200g	0.01g
	EP2102	2100g	0.01g		EX2202G	83021459	2200g	0.01g
	EP4102C	4100g	0.01g		EX4202G	83021460	4200g	0.01g
	EP4102	4100g	0.01g		EX4202G	83021460	4200g	0.01g
	EP6102C	6100g	0.01g		EX6202G	83021461	6200g	0.01g
	EP6102	6100g	0.01g		EX6202G	83021461	6200g	0.01g
	EP6101C	6100g	0.1g		EX6201G	83021462	6200g	0.1g
	EP6101	6100g	0.1g		EX6201G	83021462	6200g	0.1g
EP8101C	8100g	0.1g	EX10201G	83021464	10200g	0.1g		
EP8101	8100g	0.1g	EX10201G	83021464	10200g	0.1g		
エクスプローラーシリーズ	E1K210	12kg	0.1g	2003	EX12001G	30057115	12kg	0.1g
	E1L210	22kg	0.1g		EX24001G	30057116	22kg	0.1g
	E1M210	32kg	0.1g		EX35001G	30057117	35kg	0.1g
	E0K210	12kg	0.1g		EX12001G	30057115	12kg	0.1g
	E0L210	22kg	0.1g		EX24001G	30057116	22kg	0.1g
	E0M210	32kg	0.1g		EX35001G	30057117	35kg	0.1g
ボイジャーシリーズ	V1K210	12kg	0.1g	2003	EX12001G	30057115	12kg	0.1g
	V1L210	22kg	0.1g		EX24001G	30057116	22kg	0.1g
	V1M210	32kg	0.1g		EX35001G	30057117	35kg	0.1g
	V0K210	12kg	0.1g		EX12001G	30057115	12kg	0.1g
	V0L210	22kg	0.1g		EX24001G	30057116	22kg	0.1g
	V0M210	32kg	0.1g		EX35001G	30057117	35kg	0.1g
APシリーズ	AP110S	110g	0.1mg	1998	EX124G	83021448	120g	0.1mg
	AP210S	210g	0.1mg		EX224G	83021449	220g	0.1mg
	AP310S	310g	0.1mg		EX324G	83021450	320g	0.1mg
	AP250D	52g/210g	0.01mg/0.1mg		EX225D	30139512	120g/220g	0.1/0.01mg
Gシリーズ	G110	110g	0.1mg	1986	EX124G	83021448	120g	0.1mg
	G200	200g	0.1mg		EX224G	83021449	220g	0.1mg
	G200D	40g/200g	0.01mg/0.1mg		EX225D	30139512	120g/220g	0.1/0.01mg
ディスカバリーシリーズ	DV215CD	210/81g	0.1/0.01mg	2013	EX125D	30139510	52/120g	0.01/0.1mg
	DV114C	110g	0.1mg		EX124G	83021448	120g	0.1mg
	DV214C	210g	0.1mg		EX224G	83021449	220g	0.1mg
	DV314C	310g	0.1mg		EX324G	83021450	320g	0.1mg
エクスプローラープロシリーズ	EP12001C	12kg	0.1g	2013	EX12001G	30057115	12kg	0.1g
	EP12001	12kg	0.1g		EX12001G	30057115	12kg	0.1g
	EP22001C	22kg	0.1g		EX24001G	30057116	24kg	0.1g
	EP22001	22kg	0.1g		EX24001G	30057116	24kg	0.1g
	EP32001C	32kg	0.1g		EX35001G	30057117	35kg	0.1g
	EP32001	32kg	0.1g		EX35001G	30057117	35kg	0.1g
Eシリーズ	E10640	62g	0.1mg	2013	AX124JP	30100660	120g	0.1mg
	E11140	110g	0.1mg		AX124JP	30100660	120g	0.1mg
	E12140	210g	0.1mg		AX324JP	30100662	320g	0.1mg
	E1RR80	100g/210g	0.1mg/1mg		AX324JP	30100662	320g	0.1mg
	E12130	210g	1mg		AX223JP	30100663	220g	1mg

エクスプローラーシリーズ	E14130	410g	1mg	2003	AX423JP	30100664	420g	1mg
	E1RV70	100g/410g	1mg/0.01g		AX423JP	30100664	420g	1mg
	E16120	610g	0.01g		AX2202JP	30100667	2200g	0.01g
	E1B120	2100g	0.01g		AX2202JP	30100667	2200g	0.01g
	E1D120	4100g	0.01g		AX4202JP	30100668	4200g	0.01g
	E1RW60	1000g/4100g	0.01g/0.1g		AX4202JP	30100668	4200g	0.01g
	E1D110	4100g	0.1g		AX8201JP	30122650	8200g	0.1g
	E1F110	6100g	0.1g		AX8201JP	30122650	8200g	0.1g
	E1H110	8100g	0.1g		AX8201JP	30122650	8200g	0.1g
	E00640	62g	0.1mg		AX124JP	30100660	120g	0.1mg
	E01140	110g	0.1mg		AX124JP	30100660	120g	0.1mg
	E02140	210g	0.1mg		AX324JP	30100662	320g	0.1mg
	E0RR80	100g/210g	0.1mg/1mg		AX324JP	30100662	320g	0.1mg
	E02130	210g	1mg		AX223JP	30100663	220g	1mg
	E04130	410g	1mg		AX423JP	30100664	420g	1mg
	E0RV70	100g/410g	1mg/0.01g		AX423JP	30100664	420g	1mg
	E06120	610g	0.01g		AX1502JP	30100666	1520g	0.01g
	E0B120	2100g	0.01g		AX2202JP	30100667	2200g	0.01g
	E0D120	4100g	0.01g		AX4202JP	30100668	4200g	0.01g
	E0RW60	1000g/4100g	0.01g/0.1g		AX4202JP	30100668	4200g	0.01g
	E0D110	4100g	0.1g		AX4201JP	30122649	4200g	0.1g
	E0F110	6100g	0.1g		AX8201JP	30122650	8200g	0.1g
	E0H110	8100g	0.1g		AX8201JP	30122650	8200g	0.1g
	GTシリーズ_V1~V2	GT210	210g		1mg	1996	AX223JP	30100663
GT410		410g	1mg	AX423JP	30100664		420g	1mg
GT2100		2100g	0.01g	AX2202JP	30100667		2200g	0.01g
GT4100		4100g	0.01g	AX4202JP	30100668		4200g	0.01g
GT480		80g/400g	1mg/0.01g	AX423JP	30100664		420g	1mg
GT410D		100g/410g	1mg/0.01g	AX423JP	30100664		420g	1mg
GT4800		800g/4000g	0.01g/0.1g	AX4202JP	30100668		4200g	0.01g
GT4100D		1000g/4100g	0.01g/0.1g	AX4202JP	30100668		4200g	0.01g
GT400		400g	0.01g	STX422JP	30268871		420g	0.01g
GT4000		4100g	0.1g	AX4201JP	30122649		4200g	0.1g
GT8000		8100g	0.1g	AX8201JP	30122650		8200g	0.1g
GT8000T		8100g	0.1g	AX8201JP	30122650		8200g	0.1g
アドベンチャープロシリーズ (USBモデル)	AV114CU	110g	0.1mg	2010	AX124JP	30100660	120g	0.1mg
	AV264CU	260g	0.1mg		AX324JP	30100662	320g	0.1mg
	AV213CU	210g	1mg		AX223JP	30100663	220g	1mg
	AV413CU	410g	1mg		AX423JP	30100664	420g	1mg
	AV2102CU	2100g	0.01g		AX2202JP	30100667	2200g	0.01g
	AV4102CU	4100g	0.01g		AX4202JP	30100668	4200g	0.01g
	AV213U	210g	1mg		AX223JP	30100663	220g	1mg
	AV413U	410g	1mg		AX423JP	30100664	420g	1mg
	AV2102U	2100g	0.01g		AX2202JP	30100667	2200g	0.01g
	AV4102U	4100g	0.01g		AX4202JP	30100668	4200g	0.01g
	AV212U	210g	0.01g		STX222JP	30268870	220g	0.01g
	AV412U	410g	0.01g		STX422JP	30268871	420g	0.01g
	AV812U	810g	0.01g		AX1502JP	30100666	1520g	0.01g
	AV2101U	2100g	0.1g		AX4201JP	30122649	4200g	0.1g
	AV4101U	4100g	0.1g		AX4201JP	30122649	4200g	0.1g
	AV8101U	8100g	0.1g		AX8201JP	30122650	8200g	0.1g
AV8101CU	8100g	0.1g	AX8201JP	30122650	8200g	0.1g		
アドベンチャープロシリーズ	AV114CJP	110g	0.1mg	2013	AX124JP	30100660	120g	0.1mg
	AV264CJP	260g	0.1mg		AX324JP	30100662	320g	0.1mg
	AV213CJP	210g	1mg		AX223JP	30100663	220g	0.001g
	AV413CJP	410g	1mg		AX423JP	30100664	420g	0.001g
	AV2102CJP	2100g	0.01g		AX2202JP	30100667	2200g	0.01g
	AV4102CJP	4100g	0.01g		AX4202JP	30100668	4200g	0.01g
	AV8101CJP	8100g	0.1g		AX8201JP	30122650	8200g	0.1g
	AV213JP	210g	1mg		AX223JP	30100663	220g	0.001g
	AV413JP	410g	1mg		AX423JP	30100664	420g	0.001g
	AV2102JP	2100g	0.01g		AX2202JP	30100667	2200g	0.01g
	AV4102JP	4100g	0.01g		AX4202JP	30100668	4200g	0.01g
	AV212JP	210g	0.01g		AX622JP	30122648	620g	0.01g
	AV412JP	410g	0.01g		AX622JP	30122648	620g	0.01g
	AV812JP	810g	0.01g		AX1502JP	30100666	1500g	0.01g
	AV2101JP	2100g	0.1g		AX4201JP	30122649	4200g	0.1g
	AV4101JP	4100g	0.1g		AX4201JP	30122649	4200g	0.1g
AV8101JP	8100g	0.1g	AX8201JP	30122650	8200g	0.1g		
TSシリーズ_V1~V2	TS120S	120g	1mg	1999	PX223JP/E	30428416	220g	0.001g
	TS200S	200g	1mg		PX223JP/E	30428416	220g	0.001g
	TS400S	400g	0.01g		PX2202JP/E	30428422	2200g	0.01g
	TS400D	40g/400g	1mg/0.01g		PX2202JP/E	30428422	2200g	0.01g
	TS2KS	2000g	0.01g		PX2202JP/E	30428422	2200g	0.01g
	TS4KS	4000g	0.1g		PX4202JP/E	30428424	4200g	0.01g
	TS4KD	400g/4000g	0.01g/0.1g		PX4202JP/E	30428424	4200g	0.01g
	TPシリーズ	TP200S	200g		1mg	1996	PX223JP/E	30428416
TP400S		400g	0.01g	PX2202JP/E	30428422		2200g	0.01g
TP400D		80g/400g	1mg/0.01g	PX423JP/E	30428418		420g	0.001g
TP600S		600g	0.01g	PX2202JP/E	30428422		2200g	0.01g
TP2KS		2000g	0.01g	PX2202JP/E	30428422		2200g	0.01g
TP4KS		4000g	0.1g	PX4202JP/E	30428424		4200g	0.01g
TP4KD		800g/4000g	0.01g/0.1g	PX4202JP/E	30428424		4200g	0.01g
ASシリーズ	AS60	62g	0.1mg	1996	PX84JP/E	30428412	82g	0.1mg
	AS120	122g	0.1mg		PX124JP/E	30428413	120g	0.1mg
	AS200	202g	0.1mg		PX224JP/E	30428414	220g	0.1mg
	AS260D	62g/202g	0.1mg/1mg		PX224JP/E	30428414	220g	0.1mg
	AS64	65g	0.1mg		PX84JP	30428396	82g	0.1mg

アドベンチャーSLシリーズ	AS214	210g	0.1mg	2009	PX224JP	30428398	220g	0.1mg		
	AS153	150g	1mg		PX223JP	30428400	220g	0.001g		
	AS313	310g	1mg		PX323JP	30428401	320g	0.001g		
	AS612	610g	0.01g		PX2202JP	30428406	2200g	0.01g		
	AS1502	1500g	0.01g		PX3202JP	30428407	3200g	0.01g		
	AS3102	3100g	0.01g		PX4202JP	30428408	4200g	0.01g		
アドベンチャーシリーズ	AS6101	6100g	0.1g	2005	SPX6201JP	30268890	6200g	0.1g		
	AR0640	65g	0.1mg		PX84JP	30428396	82g	0.1mg		
	AR1140	110g	0.1mg		PX124JP	30428397	120g	0.1mg		
	AR2140	210g	0.1mg		PX224JP	30428398	220g	0.1mg		
	AR2130	210g	1mg		PX223JP	30428400	220g	0.001g		
	AR3130	310g	1mg		PX323JP	30428401	320g	0.001g		
	AR5120	510g	0.01g		PX2202JP	30428406	2200g	0.01g		
	ARB120	2100g	0.01g		PX2202JP	30428406	2200g	0.01g		
	ARC120	3100g	0.01g		PX3202JP	30428407	3200g	0.01g		
	ナビゲーターシリーズ_V1	N00330	32g		2mg	2002	NV123JP	30615798	120g	0.001g
N02120		210g	0.01g	NV222JP	30456396		220g	0.01g		
N04120		410g	0.01g	NV422JP	30456397		420g	0.01g		
N08110		810g	0.1g	NV1201JP	30456401		1200g	0.1g		
N0B110		2100g	0.1g	NV2201JP	30456402		2200g	0.1g		
N0D110		4100g	0.1g	NVT6201JP	30456405		6200g	0.1g		
N0H110		8100g	0.5g	NVT10201JP	30615804		10200g	0.1g		
N20330		32g	2mg	SPX123JP	30268880		120g	0.001g		
N22120		210g	0.01g	SPX222JP	30268882		220g	0.01g		
N24120		410g	0.01g	SPX422JP	30268883		420g	0.01g		
N28110		810g	0.1g	SPX2201JP	30268889		2200g	0.1g		
N2B110		2100g	0.1g	SPX2201JP	30268889		2200g	0.1g		
N2D110		4100g	0.1g	SPX6201JP	30268890		6200g	0.1g		
N2H110		8100g	0.5g	AX8201JP	30122650		8200g	0.1g		
PAシリーズ	PA64JP	65g	0.1mg	2018	PX84JP/E	30428412	82g	0.1mg		
	PA114JP	110g	0.1mg		PX124JP/E	30428413	120g	0.1mg		
	PA214JP	210g	0.1mg		PX224JP/E	30428414	220g	0.1mg		
	PA213JP	210g	1mg		PX223JP/E	30428416	220g	0.001g		
	PA413JP	410g	1mg		PX423JP/E	30428418	420g	0.001g		
	PA2102JP	2100g	0.01g		PX2202JP/E	30428422	2200g	0.01g		
	PA4102JP	4100g	0.01g		PX4202JP/E	30428424	4200g	0.01g		
	PA64CJP	65g	0.1mg		PX84JP	30428396	82g	0.1mg		
	PA114CJP	110g	0.1mg		PX124JP	30428397	120g	0.1mg		
	PA214CJP	210g	0.1mg		PX224JP	30428398	220g	0.1mg		
	PA213CJP	210g	1mg		PX223JP	30428400	220g	0.001g		
	PA413CJP	410g	1mg		PX423JP	30428402	420g	0.001g		
	PA2102CJP	2100g	0.01g		PX2202JP	30428406	2200g	0.01g		
	PA4102CJP	4100g	0.01g		PX4202JP	30428408	4200g	0.01g		
スカウトプロシリーズ	SP123FJP	120g	0.001g	2016	SPX123JP	30268880	120g	0.001g		
	SP2001FJP	2000g	0.1g		SPX2201JP	30268889	2200g	0.1g		
	SP202FJP	200g	0.01g		SPX222JP	30268882	220g	0.01g		
	SP4001FJP	4000g	0.1g		SPX6201JP	30268890	6200g	0.1g		
	SP6001FJP	6000g	0.1g		SPX6201JP	30268890	6200g	0.1g		
	SP401FJP	400g	0.1g		SPX421JP	30268887	420g	0.1g		
	SP402FJP	400g	0.01g		SPX422JP	30268883	420g	0.01g		
	SP6000FJP	6000g	1g		SPX8200JP	30268891	8200g	1g		
	SP601FJP	600g	0.1g		SPX621JP	30268888	620g	0.1g		
	SP602FJP	600g	0.01g		SPX622JP	30268884	620g	0.01g		
	ナビゲーターシリーズ_V2	N10330	32g		2mg	2004	NV123JP	30615798	120g	0.001g
		N12120	210g		0.01g		NV222JP	30456396	220g	0.01g
N14120		410g	0.01g	NV422JP	30456397		420g	0.01g		
N18110		810g	0.1g	NV1201JP	30456401		1200g	0.1g		
N1B110		2100g	0.1g	NV2201JP	30456402		2200g	0.1g		
N1D110		4100g	0.1g	NVT6201JP	30456405		6200g	0.1g		
N1H110		8100g	0.5g	NVT10201JP	30615804		10200g	0.1g		
N30330		32g	2mg	PX223JP	30428400		220g	0.001g		
N32120		210g	0.01g	PX2202JP	30428406		2200g	0.01g		
N34120		410g	0.01g	PX2202JP	30428406		2200g	0.01g		
N38110		810g	0.1g	PX2202JP	30428406		2200g	0.01g		
N3B110		2100g	0.1g	PX2202JP	30428406		2200g	0.01g		
N3D110		4100g	0.1g	PX4202JP	30428408		4200g	0.01g		
N3H110		8100g	0.5g	AX8201JP	30122650		8200g	0.1g		
スカウトIIシリーズ	SC2020	200g	0.01g	2003	SPX222JP	30268882	220g	0.01g		
	SC4010	400g	0.1g		SPX422JP	30268883	420g	0.01g		
	SC6010	600g	0.1g		SPX621JP	30268888	620g	0.1g		
	SCA210	1200g	0.1g		SPX2201JP	30268889	2200g	0.1g		
	SCFOA0	6000g	1g		SPX8200JP	30268891	8200g	1g		
	SR2020	200g	0.01g		SPX222JP	30268882	220g	0.01g		
	SR4010	400g	0.01g		SPX422JP	30268883	420g	0.01g		
	SR6010	600g	0.01g		SPX621JP	30268888	620g	0.1g		
	SRA210	1200g	0.1g		SPX2201JP	30268889	2200g	0.1g		
	SRFOA0	6000g	1g		SPX8200JP	30268891	8200g	1g		
スカウトシリーズ	SC2020	200g	0.01g	2000	SPX222JP	30268882	220g	0.01g		
	SC4010	400g	0.1g		SPX422JP	30268883	420g	0.01g		
	SC6010	600g	0.1g		SPX621JP	30268888	620g	0.1g		
CTシリーズ_V1~V3	CT10	10g	2mg	1998	SPX123JP	30268880	120g	0.001g		
	CT200	200g	0.01g		SPX222JP	30268882	220g	0.01g		
	CT600	600g	0.1g		SPX621JP	30268888	620g	0.1g		
	CT1200	1200g	0.1g		SPX2201JP	30268889	2200g	0.1g		
	CT6000	6000g	1g		SPX8200JP	30268891	8200g	1g		
TAシリーズ	TA103JP	100g	1mg	2013	NV123JP	30615798	120g	0.001g		
	TA152JP	150g	0.01g		NV222JP	30456396	220g	0.01g		
	TA302JP	300g	0.01g		NV422JP	30456397	420g	0.01g		

TAシリーズ	TA501JP	500g	0.1g	2019	NV622JP	30456398	620g	0.01g	
	TA1501JP	1500g	0.1g		NV2201JP	30456402	2200g	0.1g	
	TA3001JP	3000g	0.1g		NVT4201JP	30456404	4200g	0.1g	
	TA5000JP	5000g	1g		NVT6200JP	30456407	6200g	0.1g	
ナビゲーターXT(NVT)	NVT1601JP/2	1600g	0.1g	2019	NVT2201JP	30456403	2200g	0.1g	
	NVT3201JP/2	3200g	0.2g		NVT4201JP	30456404	4200g	0.1g	
	NVT6401JP/2	6400g	0.5g		NVT6201JP	30456405	6200g	0.1g	
	NVT10001JP/2	10000g	0.5g		NVT12000JP	30456408	12000g	1g	
	NVT16000JP/2	16000g	1g		NVT22000JP	30456409	22000g	1g	
	2014	NVT1601JP/3	1600g	0.1g	NVT2201JP	30456403	2200g	0.1g	
		NVT3201JP/3	3200g	0.2g	NVT4201JP	30456404	4200g	0.1g	
		NVT6401JP/3	6400g	0.5g	NVT6201JP	30456405	6200g	0.1g	
		NVT10001JP/3	10000g	0.5g	NVT12000JP	30456408	12000g	1g	
		NVT16000JP/3	16000g	1g	NVT22000JP	30456409	22000g	1g	
ナビゲーター(NV)	NV212JP/2	210g	0.01g	2019	NV222JP	30456396	220g	0.01g	
	NV511JP/2	510g	0.1g	2014	NV621JP	30456400	620g	0.1g	
	NV1101JP/2	1100g	0.1g		NV1201JP	30456401	1200g	0.1g	
	NV2101JP/2	2100g	0.1g		NV2201JP	30456402	2200g	0.1g	
	NV4101JP/2	4100g	0.2g		NVT6201JP	30456405	5100g	0.1g	
NV5101JP/2	5100g	0.5g	NVT6201JP		30456405	5100g	0.1g		
ナビゲーターXL(NVL)	NVL511JP/3	510g	0.1g	2019	NV621JP	30456400	620g	0.1g	
	NVL1101JP/3	1100g	0.1g		NV1201JP	30456401	1200g	0.1g	
	NVL2101JP/3	2100g	0.1g		NV2201JP	30456402	2200g	0.1g	
	NVL5101JP/3	5100g	0.5g		NVT6201JP	30456405	6200g	1g	
	NVL10000JP/3	10000g	1g		NVT12000JP	30456408	12000g	1g	
JE&JRシリーズ	NVL20000JP/3	20000g	1g	NVT22000JP	30456409	22000g	1g		
	JE120	120g	0.1g	2021	NV221JP	30456399	220g	0.1g	
	JE250	250g	0.1g		NV221JP	30456399	220g	0.1g	
	JE500	500g	0.1g		NV621JP	30456400	620g	0.1g	
	JR300	300g	0.1g		NV621JP	30456400	620g	0.1g	
JR2500	2500g	0.1g	NV2201JP		30456402	2200g	0.1g		
Cシリーズ	C305	300g	0.1g	2002	NV622JP	30456398	620g	0.01g	
	C505	500g	0.1g		NV622JP	30456398	620g	0.01g	
CSシリーズ	CS201JP	200g	0.1g	2019	CX221JP	30428172	220g	0.1g	
	CS2000JP	2000g	1g		CR221JP	30428177	220g	0.1g	
	CS5000JP	5000g	1g		CR621JP	30428178	620g	0.1g	
CLシリーズ	CL201JP	200g	0.1g	2013	CR2200JP	30428179	2200g	1g	
	CL5000JP	5000g	1g		CR5200JP	30428180	5200g	1g	
	CSシリーズ_V1	CS200	200g	0.1g	2005	CX221JP	30428172	220g	0.1g
		CS2000	2000g	1g		CX2200JP	30428175	2200g	1g
CS5000		5000g	2g	CX5200JP		30428176	5200g	1g	
HP&HSシリーズ	HS120	120g	0.1g	2003	PS121JP	80253216	120g	0.1g	
	HP120	120g	0.1g		PS121JP	80253216	120g	0.1g	
	HP320	320g	0.1g		PS251JP	80253217	250g	0.1g	
HHシリーズ	HH120D	60g/120g	0.1g/0.2g	2005	PS121JP	80253216	120g	0.1g	
	HH120	120g	0.1g		PS121JP	80253216	120g	0.1g	
	HH320	320g	0.1g		PS251JP	80253217	250g	0.1g	
CL-Jシリーズ	CL201J	200g	0.1g	2008	CR221JP	30428177	220g	0.1g	
	CL2000J	2000g	1g		CR2200JP	30428179	2200g	1g	
	CL5000J	5000g	2g		CR5200JP	30428180	5200g	1g	
LSシリーズ	LS200	200g	0.1g	1999	CR221JP	30428177	220g	0.1g	
	LS2000	2000g	1g		CR2200JP	30428179	2200g	1g	
	LS5000	5000g	2g		CR5200JP	30428180	5200g	1g	
	MBシリーズ	MB45	45g		0.001g/0.01%	2016	MB120	30241159	120g
MB35		35g	0.002g/0.02%	MB95	30566905		90g	0.01%/0.001g	
MB200_V1~V2	MB200	200g	0.01g	2000	MB120	30241159	120g	0.001g/0.01%	
D5000シリーズ	D51P15HR1JP	15kg	2g	2018	D52P15RTDR1JP	30479250	6/15kg	0.5/1g	
	D51P15QR1JP	15kg	2g		D52P15RTDR1JP	30479250	6/15kg	0.5/1g	
	D51P30HR1JP	30kg	2g		D52P30RTDR1JP	30479251	15/30kg	1/2g	
	D51P30QR1JP	30kg	2g		D52P30RQDL2JP	30479252	15/30kg	1/2g	
	D51P60HR1JP	60kg	5g		D52P60RTDR1JP	30479253	30/60kg	2/5g	
	D51P60HL2JP	60kg	5g		D52P60RQDL2JP	30479254	30/60kg	2/5g	
	D51P60QL2JP	60kg	5g		D52P60RQDL2JP	30479254	30/60kg	2/5g	
	D51P100HL2JP	100kg	10g		D52P150RTDL2JP	30479255	60/150kg	5/10g	
	D51P150HX2JP	150kg	10g		D52P150RTDX2JP	30479256	60/150kg	5/10g	
	D51P150QL2JP	150kg	10g		D52P150RTDL2JP	30479255	60/150kg	5/10g	
	D51P300HX2JP	300kg	20g		D52P300RTDX2JP	30479257	150/300kg	10/20g	
	D51P300QX2JP	300kg	20g		D52P300RTDX2JP	30479257	150/300kg	10/20g	
	D51XW15HR1JP	15kg	2g		2013	D52XW15WQDR6JP	30479258	6/15kg	0.5/1g
	D51XW15QR1JP	15kg	2g			D52XW15WQDR6JP	30479258	6/15kg	0.5/1g
	D51XW30HR1JP	30kg	2g			D52XW30WQDR6JP	30479259	15/30kg	1/2g
	D51XW30QR1JP	30kg	2g	D52XW30WQDR6JP		30479259	15/30kg	1/2g	
	D51XW60HL2JP	60kg	5g	D52XW60WQDL7JP		30479260	30/60kg	2/5g	
	D51XW60QL2JP	60kg	5g	D52XW60WQDL7JP		30479260	30/60kg	2/5g	
	D51XW60HR1JP	60kg	5g	D52XW60WQDL7JP		30479260	30/60kg	2/5g	
	D51XW100HL2JP	100kg	10g	D52XW150WQDX7JP		30479261	60/150kg	5/10g	
	D51XW150HX2JP	150kg	10g	D52XW150WQDX7JP		30479261	60/150kg	5/10g	
	D51XW150QL2JP	150kg	10g	D52XW150WQDX7JP		30479261	60/150kg	5/10g	
	D51XW300HX2JP	300kg	20g	D52XW300WTDX7JP		30479262	150/300kg	10/20g	
	D51XW300QX2JP	300kg	20g	D52XW300WTDX7JP		30479262	150/300kg	10/20g	
	D51XW15WR3JP	15kg	2g	2019		D52XW15WQDR6JP	30479258	6/15kg	0.5/1g
	D51XW30WR3JP	30kg	2g			D52XW30WQDR6JP	30479259	15/30kg	1/2g
	D51XW60WL4JP	60kg	5g			D52XW60WQDL7JP	30479260	30/60kg	2/5g
	D51XW150WL4JP	150kg	10g		D52XW150WQDX7JP	30479261	60/150kg	5/10g	
	D51XW300WX4JP	300kg	20g		D52XW300WTDX7JP	30479262	150/300kg	10/20g	
	CKW3R55	3kg	0.5g			30479258	6/15kg	0.5/1g	

チャンプCKWシリーズ	CKW6R55	6kg	1g	2013	D52XW15WQDR6JP	30479258	6/15kg	0.5/1g
	CKW15L55	15kg	2g		D52XW30WQDR6JP	30479259	15/30kg	1/2g
	CKW30L55	30kg	5g		D52XW30WQDR6JP	30479259	15/30kg	1/2g
チャンプSQシリーズ (CD11)	CQ10R11	10kg	1g	2009	D52P15RTDR1JP	30479250	6/15kg	0.5/1g
	CQ25R11	25kg	2g		D52P30RQDL2JP	30479252	15/30kg	1/2g
	CQ50L11	50kg	5g		D52P60RQDL2JP	30479254	30/60kg	2/5g
	CQ100L11	100kg	10g		D52P150RTDL2JP	30479255	60/150kg	5/10g
	CQ250XL11	250kg	20g		D52P300RTDX2JP	30479257	150/300kg	10/20g
チャンプSQWシリーズ (CW-11)	CQ10R11W	10kg	1g	2009	D52XW15WQDR6JP	30479258	6/15kg	0.5/1g
	CQ25R11W	25kg	2g		D52XW30WQDR6JP	30479259	15/30kg	1/2g
	CQ50L11W	50kg	5g		D52XW60WQDL7JP	30479260	30/60kg	2/5g
	CQ100L11W	100kg	10g		D52XW150WQDX7JP	30479261	60/150kg	5/10g
	CQ250XL11W	250kg	20g		D52XW300WTDX7JP	30479262	150/300kg	10/20g
D3000シリーズ	D31P30BRJP	30kg	5g	2022	D33P30B1R1JP	30683733	30kg	2g
	D31P60BRJP	60kg	10g		D33P60B1R1JP	30683734	60kg	5g
	D31P60BLJP	60kg	10g		D33P60B1L2JP	30683738	60kg	5g
	D31P150BLJP	150kg	20g		D33P150B1L2JP	30683739	150kg	10g
	D31P150BXJP	150kg	20g		D33P150B1X2JP	30683742	150kg	10g
	D31P300BXJP	300kg	50g	D33P300B1X2JP	30683743	300kg	20g	
	D32XW30BRJP	30kg	5g	2013	D33XW30C1R6JP	30683736	30kg	2g
	D32XW60BRJP	60kg	10g		D33XW60C1R6JP	30683737	60kg	5g
	D32XW60BLJP	60kg	10g		D33XW60C1L7JP	30683740	60kg	5g
	D32XW150BLJP	150kg	20g		D33XW150C1L7JP	30683741	150kg	10g
	D32XW150BXJP	150kg	20g		D33XW150C1X7JP	30683744	150kg	10g
	D32XW300BXJP	300kg	50g	D33XW300C1X7JP	30683745	300kg	20g	
	D32XW30VRJP	30kg	5g	2017	D33XW30C1R6JP	30683736	30kg	2g
	D32XW60VRJP	60kg	10g		D33XW60C1R6JP	30683737	60kg	5g
	D32XW60VLJP	60kg	10g		D33XW60C1L7JP	30683740	60kg	5g
D32XW150VLJP	150kg	20g	D33XW150C1L7JP		30683741	150kg	10g	
D32XW150VXJP	150kg	20g	D33XW150C1X7JP		30683744	150kg	10g	
D32XW300VXJP	300kg	50g	D33XW300C1X7JP	30683745	300kg	20g		
チャンプIIシリーズ	CH15R11	15kg	2g	2008	D33P15B1R1JP	30683732	15kg	1g
	CH30R11	30kg	5g		D33P30B1R1JP	30683733	30kg	2g
	CH60R11	60kg	10g		D33P60B1R1JP	30683734	60kg	5g
	CH60L11	60kg	10g		D33P60B1L2JP	30683738	60kg	5g
	CH100R11	100kg	20g		D33P150B1L2JP	30683739	150kg	10g
	CH150R11	150kg	20g		D33P150B1X2JP	30683742	150kg	10g
	CH300R11	300kg	50g		D33P300B1X2JP	30683743	300kg	20g
I-10/Lシリーズ	I-10/LB30	15kg	5g	1999	D33P30B1R1JP	30683733	30kg	2g
	I-10/LB60	30kg	10g		D33P30B1R1JP	30683733	30kg	2g
	I-10/PL100	45kg	10g		D33P60B1R1JP	30683734	60kg	5g
	I-10/PL150	70kg	20g		D33P150B1L2JP	30683739	150kg	10g
	FD3H	3kg	0.1g		V22XWE3TJP	30035490	3kg	0.5g
FDシリーズ	FD6H	6kg	0.2g	2016	V22XWE6TJP	30035491	6kg	1g
	FD15H	15kg	1g		V22XWE15TJP	30035492	15kg	2g
	FD3	3kg	0.5g		V22XWE3TJP	30035490	3kg	0.5g
	FD6	6kg	1g		V22XWE6TJP	30035491	6kg	1g
	FD15	15kg	2g		V22XWE15TJP	30035492	15kg	2g
	V21PW3JP	3kg	0.5g		V22PWE3T	30035500	3kg	0.5g
V2000シリーズ	V21PW6JP	6kg	1g	V22PWE6T	30035501	6kg	1g	
	V21PW15JP	15kg	2g	V22PWE15T	30035502	15kg	2g	
	V11P3JP	3kg	0.5g	V12P3JP	30539410	3kg	0.5g	
V1000シリーズ	V11P6JP	6kg	1g	V12P6JP	30554411	6kg	1g	
	V11P15JP	15kg	2g	V12P15JP	30554412	15kg	2g	
	V11P30JP	30kg	5g	V12P30JP	30554413	30kg	5g	
	EB3JP	3kg	0.1g	R31PE3JP	30031768	3kg	0.1g	
EBシリーズ	EB6JP	6kg	0.2g	R31PE6JP	30031769	6kg	0.2g	
	EB15JP	15kg	0.5g	R31PE15JP	30031770	15kg	0.5g	
	EB30JP	30kg	1g	R31PE30JP	30031771	30kg	1g	
	TR3RS	3kg	0.5g	R31PE3JP	30031768	3kg	0.1g	
トルーパーシリーズ	TR6RS	6kg	1g	R31PE6JP	30031769	6kg	0.2g	
	TR15RS	15kg	2g	R31PE15JP	30031770	15kg	0.5g	
	TR30RS	30kg	5g	R31PE30JP	30031771	30kg	1g	
	EC3JP	3kg	0.1g	RC31P3JP	30031808	3kg	0.1g	
	EC6JP	6kg	0.2g	RC31P6JP	30031809	6kg	0.2g	
ECシリーズ	EC15JP	15kg	0.5g	RC31P15JP	30031810	15kg	0.5g	
	EC30JP	30kg	1g	RC31P30JP	30031811	30kg	1g	
	TC3RS	3kg	0.5g	RC31P3JP	30031808	3kg	0.1g	
	TC6RS	6kg	1g	RC31P6JP	30031809	6kg	0.2g	
トルーパーカウントシリーズ	TC15RS	15kg	2g	RC31P15JP	30031810	15kg	0.5g	
	TC30RS	30kg	5g	RC31P30JP	30031811	30kg	1g	
	DS4	20kg	10g	ES50L	71138839	50kg	20g	
	DS10	50kg	20g	ES50R	71138838	50kg	20g	
DSシリーズ	DS10L	50kg	20g	ES50L	71138839	50kg	20g	
	DS20L	100kg	50g	ES100L	71138840	100kg	50g	
	DS44KL	200kg	100g	ES200L	71138823	200kg	100g	
	プリンタ_CBM910	CBM910	-	-	SF40A	30045641	-	-
	プリンタ_SF42	SF42	-	-	SF40A	30045641	-	-
C1000シリーズ	C11P9JP	9kg	5g	2023	i-C12P9JP	30820498	9kg	5g
	C11P20JP	20kg	10g		i-C12P20JP	30820499	20kg	10g
	C11P75JP	75kg	50g		i-C12P75JP	30822681	75kg	50g
SDシリーズ	SD35JP	35kg	20g	2023	i-C31M35RJP	30746011	35kg	20g
	SD75JP	75kg	50g		i-C31M75RJP	30746012	75kg	50g
	SD200JP	200kg	100g		i-C31M200RJP	30746013	200kg	100g
	SD75LJP	75kg	50g		i-C31M75LJP	30746014	75kg	50g
	SD200LJP	200kg	100g		i-C31M200LJP	30746015	200kg	100g
	ES6R	6kg	2g		i-C52X6RJP	30832025	6kg	2g
	ES30R	30kg	10g		i-C52X30RJP	30832026	30kg	10g

ESシリーズ	ES50R	50kg	20g	2023	i-C52X50RJP	30832027	50kg	20g
	ES50L	50kg	20g		i-C52X50LJP	30832028	50kg	20g
	ES100L	100kg	50g		i-C52X100LJP	30832029	100kg	50g
	ES200L	200kg	100g		i-C52X200LJP	30832030	200kg	100g
PXシリーズ	PX5202JP	5200g	0.01g	2022	PX6202JP	30694572	6200g	0.01g
	PX5202JP/E	5200g	0.01g		PX6202JP/E	30694573	6200g	0.01g